



Adult  
Ear Clip



Adult  
Finger clip



Pediatric  
Finger clip



Adult Soft  
Tip



Infant Soft  
Tip

# 1.1 Reusable SPO2 Sensor



Multi Site



Prdiatic  
Soft Tip




Adult/ Neonate/  
Ped/ Inf wrap



Neonate  
Wrap

## SpO2 Sensor & Interface Cable



### Certification

**FDA**

**CE** 0123

**NMPA**

National Medical  
Products  
Administration

### Compliance with

**EN 60601-1**

**ISO80601-2-61**

**ISO 10993-1**

**ISO 10993-5**


**ISO 10993-10**

### Specification

<b>Sensor Connector:</b>	OEM compatible, Nylon adapter with gold plated pins
<b>Cable Material &amp; Colour:</b>	TPU jacket; light gray or blue
<b>Cable Length &amp; Diameter:</b>	1.1m/1.5m/3.0m; $\phi$ 3.20mm/4.00mm or 5.00X2.50mm two line cable
<b>Sensor House :</b>	ABS shell(finger clip clip); Silicon rubber boots( soft tip and multi-site)

### Compatible with Major Patient Monitor Brands

<b>Ohmeda</b>	<b>Marquette</b>	<b>Nihon Kohden</b>	<b>Novamatrix</b>	<b>Datascope</b>	<b>Criticare</b>
<b>Philips</b>	<b>Spacelabs</b>	<b>BCI</b>	<b>Palco</b>	<b>Datex</b>	<b>Contact us for other brands</b>



	Name	Patient weight	Cable colour
<b>A</b>	Adult Finger Clip	+40kg	Gray cable
<b>B</b>	Pediatric Finger Clip	15-40kg	Gray cable
<b>C</b>	Adult Soft Tip	+40kg	Blue cable
<b>D</b>	Pediatric Soft Tip	15-40kg	Gray two line
<b>E</b>	Infant Soft Tip	3-15kg	Gray two line
<b>F</b>	Neonate Wrap	1-3kg	Gray two line
<b>G</b>	Multi-site	All patients	Gray two line
<b>H</b>	Adult ear clip	+40kg	-
<b>I</b>	Adult/Neo. Wrap	+40kg or -3kg	Gray two line
<b>J</b>	Pedia/Inf. Wrap	15-40kg	Gray two line
<b>K</b>	Multi-site	All patients	Sensormate Used with DPT tapes
<b>L</b>	Forehead/Rectal	Adult	Gray cable

## Contact Us

**EMU Medical** is here to assist in ordering the correct accessories of the Patient Monitoring System, simply email us the following information:



- Manufacturer and Model of the medical device.
- Image of the required accessories.
- Information about the part (length, connection type, etc).

We will respond promptly with all the information you need to know (part number, price, stock).

OEM Manufacture & Models	OEM REF	Connector	EMU Code	Descriptions	Ext - Cable
<b>Datascope / Mindray</b>					
Monitor Compatibilities: Passport (Datascope Tech) Accusat, Accutor 3/4 sat Passport, Accutor plus	0200-00-071-01		AC-PC-112-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-112-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-112-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-112-3	SPO2 Sensor Neonate. Wrap/3.0m	
			AC-AC-112-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-112-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-112-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-112-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-112-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-112-3	SPO2 Sensor Adult ear clip/3.0m	
	0060-00-071-010		AC-PC-120-1	SPO2 Sensor Pediatric finger clip/1.1m	Direct
			AC-PS-120-1	SPO2 Sensor Pediatric soft tip/1.1m	
			AC-IS-120-1	SPO2 Sensor Infant soft tip/1.1m	
			AC-NW-120-1	SPO2 Sensor Neonate. wrap/1.1m	
			AC-AC-120-1	SPO2 Sensor Adult finger clip/1.1m	
			AC-AS-120-1	SPO2 Sensor Adult soft tip/1.1m	
			AC-ANW-120-1	SPO2 Sensor Adu/Neo wrap/1.1m	
			AC-PIW-120-1	SPO2 Sensor Ped/Inf wrap/1.1m	
			AC-MS-120-1	SPO2 Sensor Multi-site/1.1m	
			AC-AEC-120-1	SPO2 Sensor Adult ear clip/1.1m	
Monitors: All Masimo Radical/Rainbow monitors with the latest technology. Accutor plus,EPM4, V12,V21,Accutor V, Spectrum OR Passport 2,Passport 2LT, Passport V, Passport Duo Passport Trio, Passport XG	1863/ LNCS DCI, 12653/ LNCS DBI, 2258 /LNCS YI		AC-PC-137-1	SPO2 Sensor Pediatric finger clip/1.1m	AC-EXT-114
			AC-PS-137-1	SPO2 Sensor Pediatric soft tip/1.1m	
			AC-IS-137-1	SPO2 Sensor Infant soft tip/1.1m	
			AC-NW-137-1	SPO2 Sensor Neonate. wrap/1.1m	
			AC-AC-137-1	SPO2 Sensor Adult finger clip/1.1m	
			AC-AS-137-1	SPO2 Sensor Adult soft tip/1.1m	
			AC-ANW-137-1	SPO2 Sensor Adu/Neo wrap/1.1m	
			AC-PIW-137-1	SPO2 Sensor Ped/Inf wrap/1.1m	
			AC-MS-137-1	SPO2 Sensor Multi-site/1.1m	
			AC-AEC-137-1	SPO2 Sensor Adult ear clip/1.1m	



OEM Manufacture & Models	OEM REF	Connector	EMU Code	Descriptions	Ext - Cable
<b>Datascope / Mindray</b>					
Monitor Compatibilities: Passport II/IILT, XG, Accutor Plus, Expert (Nellcor) Trio/Spectrum with Nellcor			AC-PC-126-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-126-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-126-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-126-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-126-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-126-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-126-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-126-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-126-3	SPO2 Sensor Multi-site/1.1m	
			AC-AEC-126-3	SPO2 Sensor Adult ear clip/1.1m	
<b>CSI (Criticare)</b>					
Monitor Compatibilities: Alaris IVAC: Vital Check 4400 CSI: 503DX, 503 spot, 504DX, 504Plus, 506, 506DX,507 ELC2, Scholar 508, 581,1100, 8100 Poet Plus Poet 602- 34A Scholar	934-10DN 975AD-3 940SD		AC-PC-111-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-111-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-111-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-111-3	SPO2 Sensor Neonate. Wrap/3.0m	
			AC-AC-111-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-111-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-111-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-111-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-111-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-111-3	SPO2 Sensor Adult ear clip/3.0m	
Monitor Compatibilities: 501, 503, 504, 507, 508, 581, 1100, Poet (602-3A) Scholar, 8100	934-10LN 975AD-10		AC-PC-141-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-141-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-141-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-141-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-141-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-141-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-141-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-141-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-141-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-141-3	SPO2 Sensor Adult ear clip/3.0m	



OEM Manufacture & Models	OEM REF	Connector	EMU Code	Descriptions	Ext - Cable
<b>Nonin</b>					
<p>Monitor Compatibilities: Nonin 8500, 8600, 8700, 8800 ,2500,9840 Avant 4000,9600,9700 Avant 2120PulseSense Model 7500,7500FO</p>	8000AA 8000SL		AC-PC-117-1	SPO2 Sensor Pediatric finger clip/1.1m	Direct
			AC-PS-117-1	SPO2 Sensor Pediatric soft tip/1.1m	
			AC-IS-117-1	SPO2 Sensor Infant soft tip/1.1m	
			AC-NW-117-1	SPO2 Sensor Neonate. wrap/1.1m	
			AC-AC-117-1	SPO2 Sensor Adult finger clip/1.1m	
			AC-AS-117-1	SPO2 Sensor Adult soft tip/1.1m	
			AC-ANW-117-1	SPO2 Sensor Adu/Neo wrap/1.1m	
			AC-PIW-117-1	SPO2 Sensor Ped/Inf wrap/1.1m	
			AC-MS-117-1	SPO2 Sensor Multi-site/1.1m	
			AC-AEC-117-1	SPO2 Sensor Adult ear clip/1.1m	
			AC-PC-117-3	SPO2 Sensor Pediatric finger clip/3.0m	
			AC-PS-117-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-117-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-117-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-117-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-117-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-117-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-117-3	SPO2 Sensor Ped/Inf wrap/3.0m	
AC-MS-117-3	SPO2 Sensor Multi-site/3.0m				
AC-AEC-117-3	SPO2 Sensor Adult ear clip/3.0m				
<b>GE - Critikon</b>					
<p>Monitor Compatibilities: GE Medical System: Critikon Dinamap 8700, 8710, 8720 plus, Oxyshuttle1 1 et 2</p>	9084 8731		AC-PC-118-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-118-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-118-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-118-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-118-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-118-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-118-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-118-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-118-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-118-3	SPO2 Sensor Adult ear clip/3.0m	


OEM Manufacture & Models	OEM REF	Connector	EMU Code	Descriptions	Ext - Cable
<b>GE - Critikon</b>					
Monitor Compatibilities: Storm 5300, Storm 5500 Storm 5800, Storm 5900			AC-PC-116-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-116-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-116-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-116-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-116-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-116-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-116-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-116-3	SPO2 Sensor Adult ear clip/3.0m	
			AC-MS-116-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-116-3	SPO2 Sensor Adult ear clip/3.0m	
Monitor Compatibilities: Procare, Pro Services V2, Pro 1000V3 (Nellcor Oximax Tech)			AC-PC-134-1	SPO2 Sensor Pediatric finger clip/1.1m	Direct
			AC-PS-134-1	SPO2 Sensor Pediatric soft tip/1.1m	
			AC-IS-134-1	SPO2 Sensor Infant soft tip/1.1m	
			AC-NW-134-1	SPO2 Sensor Neonate. wrap/1.1m	
			AC-AC-134-1	SPO2 Sensor Adult finger clip/1.1m	
			AC-AS-134-1	SPO2 Sensor Adult soft tip/1.1m	
			AC-ANW-134-1	SPO2 Sensor Adu/Neo wrap/1.1m	
			AC-PIW-134-1	SPO2 Sensor Ped/Inf wrap/1.1m	
			AC-MS-134-1	SPO2 Sensor Multi-site/1.1m	
			AC-AEC-134-1	SPO2 Sensor Adult ear clip/1.1m	
			AC-PC-134-3	SPO2 Sensor Pediatric finger clip/3.0m	
			AC-PS-134-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-134-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-134-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-134-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-134-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-134-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-134-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-134-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-134-3	SPO2 Sensor Adult ear clip/3.0m	



OEM Manufacture & Models	OEM REF	Connector	EMU Code	Descriptions	Ext - Cable
<b>GE Datex-Ohmeda</b>					
M-NSAT Module (AS/3,CS/3,S/5) NE-Tech Non-OxiMax & Monitor Compatibilities: Series 5 models			AC-PC-135-1	SPO2 Sensor Pediatric finger clip/1.1m	AC-EXT-151
			AC-PS-135-1	SPO2 Sensor Pediatric soft tip/1.1m	
			AC-IS-135-1	SPO2 Sensor Infant soft tip/1.1m	
			AC-NW-135-1	SPO2 Sensor Neonate. wrap/1.1m	
			AC-AC-135-1	SPO2 Sensor Adult finger clip/1.1m	
			AC-AS-135-1	SPO2 Sensor Adult soft tip/1.1m	
			AC-ANW-135-1	SPO2 Sensor Adu/Neo wrap/1.1m	
			AC-PIW-135-1	SPO2 Sensor Ped/Inf wrap/1.1m	
			AC-MS-135-1	SPO2 Sensor Multi-site/1.1m	
			AC-AEC-135-1	SPO2 Sensor Adult ear clip/1.1m	
Monitor Compatibilities: AS/3, CS/3, CH/S, CH/2S, CH-RS, CG,2G, Cardiocap 1/11/5, Capnomac Ultima, Light Monitor, M-Estpr, Oscar 11, OSP-200 Satlite Trans, Satlite Trans Plus, System 5, ULTS00, ULTS10, ULTSV1, ULTSVO	OXY-F4-N OXY-W4-N OXY-E4-N		AC-PC-113-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-113-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-113-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-113-3	SPO2 Sensor Neonate. Wrap/3.0m	
			AC-AC-113-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-113-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-113-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-113-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-113-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-113-3	SPO2 Sensor Adult ear clip/3.0m	
	OXY-F-DB		AC-PC-122-1	SPO2 Sensor Pediatric finger clip/1.1m	AC-EX-119
			AC-PS-122-1	SPO2 Sensor Pediatric soft tip/1.1m	
			AC-IS-122-1	SPO2 Sensor Infant soft tip/1.1m	
			AC-NW-122-1	SPO2 Sensor Neonate. wrap/1.1m	
			AC-AC-122-1	SPO2 Sensor Adult finger clip/1.1m	
			AC-AS-122-1	SPO2 Sensor Adult soft tip/1.1m	
			AC-ANW-122-1	SPO2 Sensor Adu/Neo wrap/1.1m	
			AC-PIW-122-1	SPO2 Sensor Ped/Inf wrap/1.1m	
			AC-MS-122-1	SPO2 Sensor Multi-site/1.1m	
			AC-AEC-122-1	SPO2 Sensor Adult ear clip/1.1m	



OEM Manufacture & Models	OEM REF	Connector	EMU Code	Descriptions	Ext - Cable
<b>GE Datex-Ohmeda</b>					
3700,3710,3770,3775, Biox3740,3760,3800 , 3900,3900P,4700, 5250RGB, Cardiocap/5, S/5,AS/3 (built in Module N-XOSAT) CS/3,(Built in module M-OSAT) (*Datex Ohmeda) M-NSAT Module (i4) NE-tech OxiMax Ohmeda Tuffsat (*Datex Ohmeda)	OXY-F1-H OXY-F4-H OXY-W4-H TS-E4-H		AC-PC-127-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-127-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-127-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-127-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-127-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-127-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-127-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-127-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-127-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-127-3	SPO2 Sensor Adult ear clip/3.0m	
			AC-PC-127-1	SPO2 Sensor Pediatric finger clip/1.1m	
			AC-PS-127-1	SPO2 Sensor Pediatric soft tip/1.1m	
			AC-IS-127-1	SPO2 Sensor Infant soft tip/1.1m	
			AC-NW-127-1	SPO2 Sensor Neonate. wrap/1.1m	
			AC-AC-127-1	SPO2 Sensor Adult finger clip/1.1m	
			AC-AS-127-1	SPO2 Sensor Adult soft tip/1.1m	
			AC-ANW-127-1	SPO2 Sensor Adu/Neo wrap/1.1m	
			AC-PIW-127-1	SPO2 Sensor Ped/Inf wrap/1.1m	
			AC-MS-127-1	SPO2 Sensor Multi-site/1.1m	
AC-AEC-127-1	SPO2 Sensor Adult ear clip/1.1m				
<b>GE - Marquette</b>					
Dash3000/4000/5000 Solar8000/8000M Solar 8000i/9500 (Tram 451M,851M) Cardiolab/Mac-Lab Corometrics 259, B20i,B40i,B60i Carescape B650, Carescape B850, Carescape V100 (*M Tech)			AC-PC-119-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-119-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-119-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-119-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-119-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-119-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-119-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-119-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-119-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-119-3	SPO2 Sensor Adult ear clip/3.0m	


OEM Manufacture & Models	OEM REF	Connector	EMU Code	Descriptions	Ext - Cable
<b>GE - Marquette</b>					
Eagle 1000/3000/4000 Dash 2000/3000/4000 Solar 8000/8000M/8000 i/9500(Tram 451/451N /851/851N).CardioLa b Mac-Lab, (Tran 451/451N) Corometrics 120/2120 Responder3000			AC-PC-138-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-138-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-138-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-138-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-138-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-138-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-138-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-138-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-138-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-138-3	SPO2 Sensor Adult ear clip/3.0m	
<b>GE - Trusignal</b>					
S/5,B20i,B40i,B60i Carescape B650, Carescape B850, Carescape V100	TS-SP-D TS-F-D TS-SA-D TS-W-D TS-SE-3 TS-W-D TS-E-D		AC-PC-125-1	SPO2 Sensor Pediatric finger clip/1.1m	AC-EX-152
			AC-PS-125-1	SPO2 Sensor Pediatric soft tip/1.1m	
			AC-IS-125-1	SPO2 Sensor Infant soft tip/1.1m	
			AC-NW-125-1	SPO2 Sensor Neonate. wrap/1.1m	
			AC-AC-125-1	SPO2 Sensor Adult finger clip/1.1m	
			AC-AS-125-1	SPO2 Sensor Adult soft tip/1.1m	
			AC-ANW-125-1	SPO2 Sensor Adu/Neo wrap/1.1m	
			AC-PIW-125-1	SPO2 Sensor Ped/Inf wrap/1.1m	
			AC-MS-125-1	SPO2 Sensor Multi-site/1.1m	
			AC-AEC-125-1	SPO2 Sensor Adult ear clip/1.1m	
	TS-F2-GE TS-SA4-GE TS-SP3-GE TS-E4-GE TS- E2-GE		AC-PC-123-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-123-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-123-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-123-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-123-3	SPO2 SensorAdult finger clip/3.0m	
			AC-AS-123-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-123-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-123-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-123-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-123-3	SPO2 Sensor Adult ear clip/3.0m	





OEM Manufacture & Models	OEM REF	Connector	EMU Code	Descriptions	Ext - Cable
<b>Mindray</b>					
Monitors: Mindray T5/T8, 1PM9800, Mindray Tech SPO2			AC-PC-115-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-115-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-115-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-115-3	SPO2 Sensor Neonate. Wrap/3.0m	
			AC-AC-115-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-115-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-115-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-115-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-115-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-115-3	SPO2 Sensor Adult ear clip/3.0m	
Monitors: Mindray T5/T8, IPM 9800, Masimo Tech SPO2 DPM 6,DPM 7 (*M tech)			AC-PC-121-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-121-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-121-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-121-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-121-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-121-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-121-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-121-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-121-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-121-3	SPO2 Sensor Adult ear clip/3.0m	
MEC-1000,MEC-1200 PM6000/7000/8000 /9000 VS 800,PM600,509B, PM9000 Express (*Mindray Tech)	512A, 512H 512H-30-79061 512G 512G-30-90607 518B 518B-30-72107 521F 512F-30-28263 512E 512E-30-90390 518C 040-000407-00 551A 0010-10-43237 ES-3212-9 0010-10-12392		AC-PC-135-1	SPO2 Sensor Pediatric finger clip/1.1m	AC-EXT-141
			AC-PS-135-1	SPO2 Sensor Pediatric soft tip/1.1m	
			AC-IS-135-1	SPO2 Sensor Infant soft tip/1.1m	
			AC-NW-135-1	SPO2 Sensor Neonate. wrap/1.1m	
			AC-AC-135-1	SPO2 Sensor Adult finger clip/1.1m	
			AC-AS-135-1	SPO2 Sensor Adult soft tip/1.1m	
			AC-ANW-135-1	SPO2 Sensor Adu/Neo wrap/1.1m	
			AC-PIW-135-1	SPO2 Sensor Ped/Inf wrap/1.1m	
			AC-MS-135-1	SPO2 Sensor Multi-site/1.1m	
			AC-F-135-1	SPO2 Sensor Forehead/1.1m	
			AC-V-135-1	SPO2 Sensor Veterinary/1.1m	
			AC-AEC-135-1	SPO2 Sensor Adult ear clip/1.1m	





OEM Manufacture & Models	OEM REF	Connector	EMU Code	Descriptions	Ext - Cable
<b>Mindray</b>					
Monitor Compatibilities: Mindray PM7000, PM6201, PM8000, PM9000 Kontron: 71388, 7840, 7845 Creative: PC9000 UP6000 Bangjian: BM9000 MEC-1000,MEC-1200 PM6000/7000/8000/9000 VS 800,PM600,509B, PM9000 Express	512A		AC-PC-124-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-124-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-124-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-124-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-124-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-124-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-124-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-124-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-124-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-124-3	SPO2 Sensor Adult ear clip/3.0m	
Monitors: Mindray T5/T8, IPM9800, VS-600/900 Nellcor Oximax Tech SPO2			AC-PC-142-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-142-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-142-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-142-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-142-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-142-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-142-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-142-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-142-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-142-3	SPO2 Sensor Adult ear clip/3.0m	
<b>CAS Medical</b>					
Monitor Compatibilities: 740-750 (Masimo Tech)	1969/ DC-12 2652/ DBI-DC8		AC-PC-116-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-116-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-116-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-116-3	SPO2 Sensor Neonate. Wrap/3.0m	
			AC-AC-116-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-116-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-116-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-116-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-116-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-116-3	SPO2 Sensor Adult ear clip/3.0m	

OEM Manufacture & Models	OEM REF	Connector	EMU Code	Descriptions	Ext - Cable
<b>NE Tech</b>					
N-10,N-20,N-30, N-85,N-100,N-180, N-190,N-200,,N-600x,NPB-190, NPB-195,NPB-3910, NPB-3920,NPB-3930, NPB-3940,NPB-4000, NPB-40,NPB-75,PM10N, PM1000N, NE-Tech Non-OxiMax	FlexMax M1192T FS-02 DS-100A FlexMax M1191T OXI-A/N+ ADH-A/N OXI-P/I+ ADH-P/ D-YS D-YSE		AC-PC-134-1	SPO2 Sensor Pediatric finger clip/1.1m	AC-EXT-188
			AC-PS-134-1	SPO2 Sensor Pediatric soft tip/1.1m	
			AC-IS-134-1	SPO2 Sensor Infant soft tip/1.1m	
			AC-NW-134-1	SPO2 Sensor Neonate. wrap/1.1m	
			AC-AC-134-1	SPO2 Sensor Adult finger clip/1.1m	
			AC-AC-134-3	SPO2 Sensor Adult finger clip/3m	
			AC-AS-134-1	SPO2 Sensor Adult soft tip/1.1m	
			AC-ANW-134-1	SPO2 Sensor Adu/Neo wrap/1.1m	
			AC-PIW-134-1	SPO2 Sensor Ped/Inf wrap/1.1m	
			AC-MS-134-1	SPO2 Sensor Multi-site/1.1m	
AC-AEC-134-1	SPO2 Sensor Adult ear clip/1.1m				
<b>Siemens</b>					
Monitor Compatibilities: Sirecust 700, Sirecust 630 Sirem SpO2 module			AC-PC-140-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-140-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-140-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-140-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-140-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-140-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-140-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-140-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-140-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-140-3	SPO2 Sensor Adult ear clip/3.0m	
			AC-PC-140-5	SPO2 Sensor Pediatric finger clip/5.0m	
			AC-PS-140-5	SPO2 Sensor Pediatric soft tip/5.0m	
			AC-IS-140-5	SPO2 Sensor Infant soft tip/5.0m	
			AC-NW-140-5	SPO2 Sensor Neonate. wrap/5.0m	
			AC-AC-140-5	SPO2 Sensor Adult finger clip/5.0m	
			AC-AS-140-5	SPO2 Sensor Adult soft tip/5.0m	
			AC-ANW-140-5	SPO2 Sensor Adu/Neo wrap/5.0m	
			AC-PIW-140-5	SPO2 Sensor Ped/Inf wrap/5.0m	
AC-MS-140-5	SPO2 Sensor Multi-site/5.0m				
AC-AEC-140-5	SPO2 Sensor Adult ear clip/5.0m				

OEM Manufacture & Models	OEM REF	Connector	EMU Code	Descriptions	Ext - Cable
<b>Siemens</b>					
Monitor Compatibilities: SC6002XL, 7000, SC8000, SC9000XL			AC-PC-147-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct ****
			AC-PS-147-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-147-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-147-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-147-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-147-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-147-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-147-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-147-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-147-3	SPO2 Sensor Adult ear clip/3.0m	
			AC-PC-147-5	SPO2 Sensor Pediatric finger clip/5.0m	
			AC-PS-147-5	SPO2 Sensor Pediatric soft tip/5.0m	
			AC-IS-147-5	SPO2 Sensor Infant soft tip/5.0m	
			AC-NW-147-5	SPO2 Sensor Neonate. wrap/5.0m	
			AC-AC-147-5	SPO2 Sensor Adult finger clip/5.0m	
			AC-AS-147-5	SPO2 Sensor Adult soft tip/5.0m	
			AC-ANW-147-5	SPO2 Sensor Adu/Neo wrap/5.0m	
			AC-PIW-147-5	SPO2 Sensor Ped/Inf wrap/5.0m	
			AC-MS-147-5	SPO2 Sensor Multi-site/5.0m	
AC-AEC-147-5	SPO2 Sensor Adult ear clip/5.0m				
<b>Masimo (M Tech)</b>					
	1864/ LNCS® DCIP 1863/ LNCS® DCI 2653/ LNCS® DBI 2258/ LNCS® YI 1895/ LNCS® TC-I		AC-PC-137-1	SPO2 Sensor Pediatric finger clip/1.1m	AC-EXT-115 AC-EXT-159
			AC-PS-137-1	SPO2 Sensor Pediatric soft tip/1.1m	
			AC-IS-137-1	SPO2 Sensor Infant soft tip/1.1m	
			AC-NW-137-1	SPO2 Sensor Neonate. wrap/1.1m	
			AC-AC-137-1	SPO2 Sensor Adult finger clip/1.1m	
			AC-AS-137-1	SPO2 Sensor Adult soft tip/1.1m	
			AC-ANW-137-1	SPO2 Sensor Adu/Neo wrap/1.1m	
			AC-PIW-137-1	SPO2 Sensor Ped/Inf wrap/1.1m	
			AC-MS-137-1	SPO2 Sensor Multi-site/1.1m	
			AC-AEC-137-1	SPO2 Sensor Adult ear clip/1.1m	




OEM Manufacture & Models	OEM REF	Connector	EMU Code	Descriptions	Ext - Cable
<b>Masimo (M Tech)</b>					
Compatible With M Tech Module 6pin connector	1276/ LNOP® DCIP 1269/ LNOP® DCI 1544/ LNOP® YI 1794/ LNOP® TC-I		AC-PC-128-1	SPO2 Sensor Pediatric finger clip/1.1m	AC-EXT-115 AC-EXT-159
			AC-PS-128-1	SPO2 Sensor Pediatric soft tip/1.1m	
			AC-IS-128-1	SPO2 Sensor Infant soft tip/1.1m	
			AC-NW-128-1	SPO2 Sensor Neonate. wrap/1.1m	
			AC-AC-128-1	SPO2 Sensor Adult finger clip/1.1m	
			AC-AS-128-1	SPO2 Sensor Adult soft tip/1.1m	
			AC-ANW-128-1	SPO2 Sensor Adu/Neo wrap/1.1m	
			AC-PIW-128-1	SPO2 Sensor Ped/Inf wrap/1.1m	
			AC-MS-128-1	SPO2 Sensor Multi-site/1.1m	
			AC-AEC-128-1	SPO2 Sensor Adult ear clip/1.1m	
			AC-PC-128-3	SPO2 Sensor Pediatric finger clip/3.0m	
			AC-PS-128-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-128-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-128-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-128-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-128-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-128-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-128-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-128-3	SPO2 Sensor Multi-site/3.0m	
AC-AEC-128-3	SPO2 Sensor Adult ear clip/3.0m				
Compatible With M Tech Module 14pin connector			AC-PC-116-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-116-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-116-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-116-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-116-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-116-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-116-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-116-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-116-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-116-3	SPO2 Sensor Adult ear clip/3.0m	

OEM Manufacture & Models	OEM REF	Connector	EMU Code	Descriptions	Ext - Cable
<b>Masimo (M Tech)</b>					
Monitors: All Masimo Radical/Rainbow monitors with the latest technology.Has chips inside!			AC-PC-133-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-133-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-133-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-133-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-133-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-133-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-133-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-133-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-133-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-133-3	SPO2 Sensor Adult ear clip/3.0m	
1396/1969/ DC8 Masimo Metal Connector			AC-PC-116-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-116-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-116-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-116-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-116-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-116-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-116-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-116-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-116-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-116-3	SPO2 Sensor Adult ear clip/3.0m	
<b>Physio - Control - Masimo</b>					
Monitor Compatibilities: Lifepak 12, 15, 20 etc, Masimo Tech.			AC-PC-116-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-116-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-116-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-116-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-116-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-116-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-116-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-116-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-116-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-116-3	SPO2 Sensor Adult ear clip/3.0m	

OEM Manufacture & Models	OEM REF	Connector	EMU Code	Descriptions	Ext - Cable
<b>Physio - Control - Masimo</b>					
Monitor Compatibilities: Lifepak 12			AC-PC-149-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-149-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-149-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-149-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-149-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-149-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-149-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-149-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-149-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-149-3	SPO2 Sensor Adult ear clip/3.0m	
<b>Philips</b>					
Monitor Compatibilities:783 52A/C, 78354A/C, 78834C, M1020A, M1025A/B, M1025A, M1350B, M1722A/B, M1723A/B, M2475B ,M1020A	M1193T M1194A M1190A		AC-PC-130-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct ***
			AC-PS-130-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-130-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-130-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-130-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-130-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-130-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-130-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-130-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-130-3	SPO2 Sensor Adult ear clip/3.0m	
			AC-PC-130-5	SPO2 Sensor Pediatric finger clip/5.0m	
			AC-PS-130-5	SPO2 Sensor Pediatric soft tip/5.0m	
			AC-IS-130-5	SPO2 Sensor Infant soft tip/5.0m	
			AC-NW-130-5	SPO2 Sensor Neonate. wrap/5.0m	
			AC-AC-130-5	SPO2 Sensor Adult finger clip/5.0m	
			AC-AS-130-5	SPO2 Sensor Adult soft tip/5.0m	
			AC-ANW-130-5	SPO2 Sensor Adu/Neo wrap/5.0m	
			AC-PIW-130-5	SPO2 Sensor Ped/Inf wrap/5.0m	
			AC-MS-130-5	SPO2 Sensor Multi-site/5.0m	
			AC-AEC-130-5	SPO2 Sensor Adult ear clip/5.0m	




OEM Manufacture & Models	OEM REF	Connector	EMU Code	Descriptions	Ext - Cable
<b>Philips</b>					
Monitor Compatibilities: MP20/30/40/ PM50/60/70/90 VN4,VM6,VM8(*M tech)			AC-PC-114-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-114-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-114-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-114-3	SPO2 Sensor Neonate. Wrap/3.0m	
			AC-AC-114-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-114-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-114-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-114-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-114-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-114-3	SPO2 Sensor Adult ear clip/3.0m	
Monitor Compatibilities:783 52B, 78354B, 78352A,78352C,788 34C, 78354A, 78354C,78834C, M1020A, ,M1025B,M1025A, M1350B/C, M1722A/B, M1723A/B, M2475B, M2601A, M300A, M3001A, M3500B, M4735A	M1192A M1192T M1194A M1193A M1195A M1196A M1191A M1193T		AC-PC-146-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct ****
			AC-PS-146-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-146-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-146-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-146-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-146-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-146-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-146-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-146-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-146-3	SPO2 Sensor Adult ear clip/3.0m	
			AC-PC-146-5	SPO2 Sensor Pediatric finger clip/5.0m	
			AC-PS-146-5	SPO2 Sensor Pediatric soft tip/5.0m	
			AC-IS-146-5	SPO2 Sensor Infant soft tip/5.0m	
			AC-NW-146-5	SPO2 Sensor Neonate. wrap/5.0m	
			AC-AC-146-5	SPO2 Sensor Adult finger clip/5.0m	
			AC-AS-146-5	SPO2 Sensor Adult soft tip/5.0m	
			AC-ANW-146-5	SPO2 Sensor Adu/Neo wrap/5.0m	
			AC-PIW-146-5	SPO2 Sensor Ped/Inf wrap/5.0m	
			AC-MS-146-5	SPO2 Sensor Multi-site/5.0m	
			AC-AEC-146-5	SPO2 Sensor Adult ear clip/5.0m	

OEM Manufacture & Models	OEM REF	Connector	EMU Code	Descriptions	Ext - Cable	
<b>Coreplus</b>						
			AC-PC-134-1	SPO2 Sensor Pediatric finger clip/1.1m	Direct	
			AC-PS-134-1	SPO2 Sensor Pediatric soft tip/1.1m		
			AC-IS-134-1	SPO2 Sensor Infant soft tip/1.1m		
			AC-NW-134-1	SPO2 Sensor Neonate. wrap/1.1m		
			AC-AC-134-1	SPO2 Sensor Adult finger clip/1.1m		
			AC-AS-134-1	SPO2 Sensor Adult soft tip/1.1m		
			AC-ANW-134-1	SPO2 Sensor Adu/Neo wrap/1.1m		
			AC-PIW-134-1	SPO2 Sensor Ped/Inf wrap/1.1m		
			AC-MS-134-1	SPO2 Sensor Multi-site/1.1m		
			AC-AEC-134-1	SPO2 Sensor Adult ear clip/1.1m		
			AC-PC-134-3	SPO2 Sensor Pediatric finger clip/3.0m		
			AC-PS-134-3	SPO2 Sensor Pediatric soft tip/3.0m		
			AC-IS-134-3	SPO2 Sensor Infant soft tip/3.0m		
			AC-NW-134-3	SPO2 Sensor Neonate. wrap/3.0m		
			AC-AC-134-3	SPO2 Sensor Adult finger clip/3.0m		
			AC-AS-134-3	SPO2 Sensor Adult soft tip/3.0m		
			AC-ANW-134-3	SPO2 Sensor Adu/Neo wrap/3.0m		
			AC-PIW-134-3	SPO2 Sensor Ped/Inf wrap/3.0m		
			AC-MS-134-3	SPO2 Sensor Multi-site/3.0m		
			AC-AEC-134-3	SPO2 Sensor Adult ear clip/3.0m		
<b>NIHON KOHDEN</b>						
Lifescope S BSS-9800 Lifescope M BSM 9510 Lifescope P BSM 4100 Lifescope P BSM 4101 Lifescope P BSM 4103/4111 Lifescope P BSM 4113 Life scope A BSM 5100 Life scope A BSM 5105 Life scope A BSM 5135 Life scopeTR BSM 6000			AC-PC-132-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct	
				AC-PS-132-3		SPO2 Sensor Pediatric soft tip/3.0m
				AC-IS-132-3		SPO2 Sensor Infant soft tip/3.0m
				AC-NW-132-3		SPO2 Sensor Neonate. wrap/3.0m
				AC-AC-132-3		SPO2 Sensor Adult finger clip/3.0m
				AC-AS-132-3		SPO2 Sensor Adult soft tip/3.0m
				AC-ANW-132-3		SPO2 Sensor Adu/Neo wrap/3.0m
				AC-PIW-132-3		SPO2 Sensor Ped/Inf wrap/3.0m
				AC-MS-132-3		SPO2 Sensor Multi-site/3.0m
				AC-AEC-132-3		SPO2 Sensor Adult ear clip/3.0m

OEM Manufacture & Models	OEM REF	Connector	EMU Code	Descriptions	Ext - Cable
<b>NIHON KOHDEN</b>					
Life scopeTR, BSM 6301K Life scopeTR, BSM 6501K Life scopeTR, BSM 6701K Lifescope L BSM 2350 LifescopeL BSM 2351 LifescopeL BSM 2353 Lifescopei BSM 2300 Lifescopei BSM 2301A Lifescopei BSM 2301K, Lifescopei BSM 2303K LifescopeL BSM 2351 VISMO PVM-2701/PVM-2703	TL-201T A/SPO2-ADVT TL-220T		AC-PC-129-1	SPO2 Sensor Pediatric finger clip/1.5m	AC-EXT-138
			AC-PS-129-1	SPO2 Sensor Pediatric soft tip/1.5m	
			AC-IS-129-1	SPO2 Sensor Infant soft tip/1.5m	
			AC-NW-129-1	SPO2 Sensor Neonate. wrap/1.5m	
			AC-AC-129-1	SPO2 Sensor Adult finger clip/1.5m	
			AC-AS-129-1	SPO2 Sensor Adult soft tip/1.5m	
			AC-ANW-129-1	SPO2 Sensor Adu/Neo wrap/1.5m	
			AC-PIW-129-1	SPO2 Sensor Ped/Inf wrap/1.5m	
			AC-MS-129-1	SPO2 Sensor Multi-site/1.5m	
			AC-AEC-129-1	SPO2 Sensor Adult ear clip/1.5m	
OPV-1500			AC-PC-136-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-136-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-136-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-136-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-136-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-136-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-136-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-136-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-136-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-136-3	SPO2 Sensor Adult ear clip/3.0m	
BSM-1101, BSM-2101 BSM-2301, BSM-7100, BSM-7200, OLV-1100A			AC-PC-139-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-139-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-139-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-139-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-139-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-139-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-139-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-139-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-139-3	SPO2 Sensor Multi-site/3.0m	
			AC-AEC-139-3	SPO2 Sensor Adult ear clip/3.0m	



OEM Manufacture & Models	OEM REF	Connector	EMU Code	Descriptions	Ext - Cable
<b>Novamatrix</b>					
<p>Monitor            Compatibilities:            Novamatrix: 505, 510, 511, 515, 515A, 515B, 515C, 520, 520A, 7000, 7100, Oxypleth;            Spacelabs: 90351-D/G, 90369, 9046, 90466, 90467, 90489, 90651-07/08, 90419-01, 90651A-07; Nico Cardiopulmonary.</p>			AC-PC-131-3	SPO2 Sensor Pediatric finger clip/3.0m	Direct
			AC-PS-131-3	SPO2 Sensor Pediatric soft tip/3.0m	
			AC-IS-131-3	SPO2 Sensor Infant soft tip/3.0m	
			AC-NW-131-3	SPO2 Sensor Neonate. wrap/3.0m	
			AC-AC-131-3	SPO2 Sensor Adult finger clip/3.0m	
			AC-AS-131-3	SPO2 Sensor Adult soft tip/3.0m	
			AC-ANW-131-3	SPO2 Sensor Adu/Neo wrap/3.0m	
			AC-PIW-131-3	SPO2 Sensor Ped/Inf wrap/3.0m	
			AC-MS-131-3	SPO2 Sensor Adult ear clip/3.0m	
			AC-AEC-131-3	SPO2 Sensor Multi-site/3.0m	
			AC-PC-131-1	SPO2 Sensor Pediatric finger clip/1.1m	
			AC-PS-131-1	SPO2 Sensor Pediatric soft tip/1.1m	
			AC-IS-131-1	SPO2 Sensor Infant soft tip/1.1m	
			AC-NW-131-1	SPO2 Sensor Neonate. wrap/1.1m	
			AC-AC-131-1	SPO2 Sensor Adult finger clip/1.1m	
			AC-AS-131-1	SPO2 Sensor Adult soft tip/1.1m	
			AC-ANW-131-1	SPO2 Sensor Adu/Neo wrap/1.1m	
			AC-PIW-131-1	SPO2 Sensor Ped/Inf wrap/1.1m	
			AC-MS-131-1	SPO2 Sensor Multi-site/1.1m	
			AC-AEC-131-1	SPO2 Sensor Adult ear clip/1.1m	

## 1.2 SPO2 Interface Cable



OEM REF	EMU Ref	Descriptions	Monitor end Connector	Sensor end Connector
<b>CSI/Criticare</b>				
518DD & 518LD	AC-EXT-110	SPO2 adapter cable 2.2m		
	AC-EXT-111	SPO2 adapter cable 2.2m		
Monitor Compatibilities: 506 Spot, 507SD, 503 Spot				
<b>Datascope</b>				
0012-00-0516-02	AC-EXT-112	SPO2 adapter cable 2.2m		
Monitor Compatibilities: Datascope: Accusat, Accutorr 3 SAT, Accutorr 4 SAT, Accutorr Plus, Passport, Passport 2, Passport 2LT, Passport IL, Passport XL and Passport XG				
<b>Dolphin</b>				
110	AC-EXT-113	SPO2 adapter cable 2.2m		
Monitor Compatibilities: Dolphin ONE oximeters				
<b>Masimo(M Tech)</b>				
LNCS Patient Cable, LNC-4(REF 2017), LNC-10(REF 1814), LNC-14(REF 2013)	AC-EXT-114	SPO2 adapter cable M Tech to LNCS 2.2m		
PC 04, PC 08, PC 12	AC-EXT-115	SPO2 adapter cable M Tech LNOP non Radical 3.0m		
	AC-EXT-159	SPO2 adapter cable M Tech LNOP non Radical 2.2m		
<b>Drager/Siemens</b>				
4533840	AC-EXT-116	SPO2 adapter cable Non Oximax 2.2m		
3375834	AC-EXT-117	SPO2 adapter cable Non Oximax 2.2m		
Monitor Compatibilities: Oxipac 2000, Oxipac 2500				
<b>Nonin</b>				
E708-080	AC-EXT-118	SPO2 adapter cable 2.2m		
Monitor Compatibilities: Nonin 8604, 8604D				











OEM REF	EMU Ref	Descriptions	Monitor end Connector	Sensor end Connector
<b>GE-Datex-Ohmeda</b>				
OXY-C3	AC-EXT-119	SPO2 adapter cable Datex 2.2m		
Monitor Compatibilities: S/5, AS/3, CS/3 Modular Monitors with M-series modules, Cardiocap/5, Cardiocap S/5 Light Monitor, Capnomac Ultima, Satellite Plus II				
2006644-001 E9004GE	AC-EXT-120	SPO2 adapter cable NE-Tech Non-OxiMax 2.2m		
	AC-EXT-121	SPO2 adapter cable NE-Tech Non-OxiMax 3.0m		
Monitor Compatibilities: Marquette Dash 2000, Dash 3000, Eagle Series, SOLAR Series, 7000 Series, Tram Modular Series				
LNC-10-GE (REF 2016) 2027263-002	AC-EXT-122	SPO2 adapter cable Masimo Tech 2.2m		
	AC-EXT-124	SPO2 adapter cable Masimo Tech 3.0m		
Monitor Compatibilities: DASH, SOLAR Series (Masimo Model)				
2021406-001	AC-EXT-137	SPO2 adapter cable Nellcor Oximax 3.0m		
Monitor Compatibilities: Oximax Model (for Nellcor Oximax probe)				
	AC-EXT-123	SPO2 adapter cable for use with disposable Masimo LNOP Sensors 2.2m		
TS-G3	AC-EXT-152	SpO2 Adapter Cable GE Trusignal Sensor 2.2m 0.3m		
TS-H3	AC-EXT-125	SPO2 adapter cable GE Trusignal Sensor 3.0m		
TS-N3	AC-EXT-126	SpO2 Adapter Cable GE Trusignal Sensor 2.2m		
TS-M3	AC-EXT-127	SpO2 Adapter Cable GE Trusignal Sensor 2.2m		


OEM REF	EMU Ref	Descriptions	Monitor end Connector	Sensor end Connector
<b>Philips</b>				
M1900A/B	AC-EXT-128	SPO2 adapter cable M tech LNOP 2.2m		
Monitor Compatibilities: 78352C, 78354C, 78834C, M1020A, M1025B, M1205A, M1350, M1722A/B, M1732A/B, M2475B, Viridia A1, Viridia A3				
M1943A	AC-EXT-129	SPO2 adapter cable NE-Tech Non-OxiMax 2.2m		
M1943AL	AC-EXT-130	SPO2 adapter cable NE-Tech Non-OxiMax 3.0m		
Monitor Compatibilities: Philips Medical: Heartstream XL, Intelliview MP70, M2601A, M3000A, M3500B, M4735A, M1020B, M3000A/M3001A MMS, M3002A MMS X2, M8102A MP2, M8105A MP5				
M1940A/B	AC-EXT-131	SPO2 adapter cable 2.2m		
Monitor Compatibilities: M1020A, M2475B, M1732A/B, M1722A/B, M1020A, 78833C, 78833A, 78834A, 78352C, 78352A5				
M1943NL	AC-EXT-135	SPO2 adapter cable NE-Tech Non-OxiMax 3.0m		
Monitor Compatibilities: All new Philips Oximax Monitors eg mms X2, MX800 etc.				
M1941A/B	AC-EXT-132	SPO2 adapter cable (Philips SPO2 cable)2.2m		
LNOP® PC-12 P8 REF 1847	AC-EXT-133	SPO2 adapter cable Masimo Tech - LNOP2.2m		
LNC-10 P8 REF 1902	AC-EXT-134	SPO2 adapter cable Masimo LNCS 2.2m		
<b>EDAN</b>				
M13-30691 (M9,M9A,M8,M8A, M8B,M3,M3B)	AC-EXT-141	SPO2 adapter cable EDAN Oximetry 2.2m		
	AC-EXT-142	SPO2 adapter cable EDAN Oximetry 3.0m		
Monitor Compatibilities: Mindray PM-600, MEC 100, MEC 200, PM 7000, 6201, PM 8000, PM 9000, VS-800, Datascope Trio, Accutorr V, Passport V				

OEM REF	EMU Ref	Descriptions	Monitor end Connector	Sensor end Connector
<b>Nihon Kohden</b>				
JL-650P	AC-EXT-136	SpO2 Adapter Cable NE-Tech OxiMax 3.0m		
JL-631P	AC-EXT-139	SpO2 Adapter Cable M Tech 3.0m		
JL-900P BSM/PVM	AC-EXT-138	SpO2 Adapter Cable NK 3.0m		
JL-302T OPV- 1500,OLV	AC-EXT-140	SpO2 Adapter Cable NK 3.0m		
<b>Mindray</b>				
0010-20-42594, Accutorr V, DPM 3, DPM 4, DPM5, Passport V	AC-EXT-141	SPO2 adapter cable NE-Tech Non-OxiMax 2.2m		
	AC-EXT-142	SPO2 adapter cable NE-Tech Non-OxiMax 3.0m		
0010-30-42737 0010-20-42710 0515-30- 11221	AC-EXT-143	SPO2 adapter cable Non Oximax 2.2m		
M technology 0010-30-12452	AC-EXT-144	SPO2 adapter cable NE-Tech Non-OxiMax 2.2m		
	AC-EXT-145	SPO2 adapter cable Non Oximax 2.2m		
(0012-00-1652),1.2m (0012-00-1599),3.1m (0012-00-1653),4.3m Accutorr V, Duo, Trio, Passport 2, PassportV, Spectrum/OR, V12, V21				
0010-20-42710 T5/T8,DPM 6, DPM 7	AC-EXT-146	SPO2 adapter cable NE-Tech Non-OxiMax 2.2m		
<b>Novamatrix</b>				
8776-00	AC-EXT-154	SPO2 adapter cable 2.2m		
Monitor Compatibilities: Novamatrix 505, 510, 511, 515B/C 520, 7000, 7100				



OEM REF	EMU Ref	Descriptions	Monitor end Connector	Sensor end Connector
<b>Nellcor Oximax (NE Tech)</b>				
DOC-10 + Oximax adaptor	AC-EXT-147	N550 Oximax Technology SPO2 adapter cable 3.0m		
Monitor Compatibilities: N-595 and N-600; CAS Apnea Monitor; Cardell Veterinary, 9500HD; Colin BP-S510 and T105NXM; Digicare LW600; Dinamap Carescape V100, ProCare, Pro Series, ProSeries V20, Pro 300, Pro1000, Pro 1000 VS, E-NSAT Modules for i4 Modular Monitors with Oximax option; Fukuda Denshi DS-7100; GE Procure; Omron; Welch Allyn Spot Vital Signs 42NTB-E1 LXI.				
39783	AC-EXT-148	Oximax Technology SPO2 adapter cable NE-Tech Non-OxiMax 0.3m		
	AC-EXT-149	Oximax Technology SPO2 adapter cable NE-Tech Non-OxiMax 2.2m		
	AC-EXT-150	Oximax Technology SPO2 adapter cable NE-Tech Non-OxiMax 3.0m		
Monitor Compatibilities: Nellcor suitable for Oximax technology SPO2 sensor, Welch Allyn				
SCP-10-MC-10	AC-EXT-151	SPO2 adapter cable Oxismart NE-Tech Non-OxiMax 2.2m		
Monitor Compatibilities: Colin: BP88/306; Critikon Dinamap ProCare, Pro Series, Pro 300, Pro 1000, M-NSAT Modules; Datascope: Accutorr Plus, Duo, Passport 2, Passport XG, Spectrum, Trio; Nellcor: NPB-290, NPB-295, N-390, N-395, NPB-290, NPB-295, N-390, N-395, N-3000, N-6000, Symphony N-3000; Nihon Kohden: BSM-2300, BSM-4100, BSM-5100, BSM-4104A, BSM-4114A; Ohmeda enGuard CM5; Physio Control Lifepak 12.				
<b>BCI</b>				
3311	AC-EXT-155	SPO2 adapter cable		
Monitor Compatibilities: 6100, 9100 Advisor, Autocorr, Mini Torr, Mini Torr Plus				

# 1.3 Disposable Spo2 Sensor









OEM Part Number	Textile Adhesive REF#	Foam Adhesive REF#	Non-Adhesive REF#	Sensor EndConnector	Patient Size
<b>General</b>					
	CO-T-A42	CO-F-A42	CO-N-A42		Adult(>30kg)
	CO-T-P42	CO-F-P42	CO-N-P42		Pediatric(15-40kg)
	CO-T-I42	CO-F-I42	CO-N-I42		Infant(1-15kg)
	CO-T-N42	CO-F-N42	CO-N-N42		Neonate(<3kg)
<b>NE tech OxiSmart</b>					
D25	CO-T-A01	CO-F-A01	CO-N-A01		Adult(>30kg)
D20	CO-T-P01	CO-F-P01	CO-N-P01		Pediatric(15-40kg)
I20	CO-T-I01	CO-F-I01	CO-N-I01		Infant(1-15kg)
N25	CO-T-N01	CO-F-N01	CO-N-N01		Neonate(<3kg)
<b>NE tech OxiMax</b>					
MAX-A	CO-T-A10	CO-F-A10	CO-N-A10		Adult(>30kg)
MAX-P	CO-T-P10	CO-F-P10	CO-N-P10		Pediatric(15-40kg)
MAX-I	CO-T-I10	CO-F-I10	CO-N-I10		Infant(1-15kg)
MAX-N	CO-T-N10	CO-F-N10	CO-N-N10		Neonate(<3kg)
<b>Novamatrix</b>					
AS110	CO-T-A03	CO-F-A03	CO-N-A03		Adult(>30kg)
AS120	CO-T-P03	CO-F-P03	CO-N-P03		Pediatric(15-40kg)
AS120	CO-T-I03	CO-F-I03	CO-N-I03		Infant(1-15kg)
AS140	CO-T-N03	CO-F-N03	CO-N-N03		Neonate(<3kg)
<b>Criticare</b>					
570SD	CO-T-A05	CO-F-A05	CO-N-A05		Adult(>30kg)
571SD	CO-T-P05	CO-F-P05	CO-N-P05		Pediatric(15-40kg)
572SD	CO-T-I05	CO-F-I05	CO-N-I05		Infant(1-15kg)
573SD	CO-T-N05	CO-F-N05	CO-N-N05		Neonate(<3kg)
<b>BCI</b>					
1300	CO-T-A06	CO-F-A06	CO-N-A06		Adult(>30kg)
1301	CO-T-P06	CO-F-P06	CO-N-P06		Pediatric(15-40kg)
1303	CO-T-I06	CO-F-I06	CO-N-I06		Infant(1-15kg)
1302	CO-T-N06	CO-F-N06	CO-N-N06		Neonate(<3kg)
<b>Datascope</b>					
0998-00-0076-05	CO-T-A07	CO-F-A07	CO-N-A07		Adult(>30kg)
0998-00-0076-04	CO-T-P07	CO-F-P07	CO-N-P07		Pediatric(15-40kg)
0998-00-0074-03	CO-T-I07	CO-F-I07	CO-N-I07		Infant(1-15kg)
0998-00-0074-04	CO-T-N07	CO-F-N07	CO-N-N07		Neonate(<3kg)
<b>Nonin</b>					
6000CA/7000A	CO-T-A08	CO-F-A08	CO-N-A08		Adult(>30kg)
6000CP/7000P	CO-T-P08	CO-F-P08	CO-N-P08		Pediatric(15-40kg)
6000CI/7000I	CO-T-I08	CO-F-I08	CO-N-I08		Infant(1-15kg)
6000CN/7000N	CO-T-N08	CO-F-N08	CO-N-N08		Neonate(<3kg)



OEM Part Number	Textile Adhesive REF#	Foam Adhesive REF#	Non-Adhesive REF#	Sensor EndConnector	Patient Size
<b>GE Ohmeda</b>					
SAS-AF	CO-T-A02	CO-F-A02	CO-N-A02		Adult(>30kg)
SAS-AP	CO-T-P02	CO-F-P02	CO-N-P02		Pediatric(15-40kg)
SAS-AP	CO-T-I02	CO-F-I02	CO-N-I02		Infant(1-15kg)
SAS-AF	CO-T-N02	CO-F-N02	CO-N-N02		Neonate(<3kg)
<b>GE Datex</b>					
SAS-AF	CO-T-A09	CO-F-A09	CO-N-A09		Adult(>30kg)
SAS-AP	CO-T-P09	CO-F-P09	CO-N-P09		Pediatric(15-40kg)
SAS-AP	CO-T-I09	CO-F-I09	CO-N-I09		Infant(1-15kg)
SAS-AF	CO-T-N09	CO-F-N09	CO-N-N09		Neonate(<3kg)
<b>GE TruSignal+</b>					
TS-AF-10	CO-T-A26	CO-F-A26	CO-N-A26		Adult(>30kg)
TS-AF-10	CO-T-P26	CO-F-P26	CO-N-P26		Pediatric(15-40kg)
TS-AF-10	CO-T-I26	CO-F-I26	CO-N-I26		Infant(1-15kg)
TS-AF-10	CO-T-N26	CO-F-N26	CO-N-N26		Neonate(<3kg)
<b>GE OxyTip+</b>					
OXY-AF-10	CO-T-A127	CO-F-A127	CO-N-A127		Adult(>30kg)
OXY-AP-10	CO-T-P127	CO-F-P127	CO-N-P127		Pediatric(15-40kg)
OXY-AF-10	CO-T-I127	CO-F-I127	CO-N-I127		Infant(1-15kg)
OXY-AF-10	CO-T-N127	CO-F-N127	CO-N-N127		Neonate(<3kg)
<b>GE M tech Rainbow</b>					
LNCS Amtx	CO-T-A49	CO-F-A49	CO-N-A49		Adult(>30kg)
LNCS Pmtx	CO-T-P49	CO-F-P49	CO-N-P49		Pediatric(15-40kg)
LNCS Inf	CO-T-I49	CO-F-I49	CO-N-I49		Infant(1-15kg)
LNCS Neo	CO-T-N49	CO-F-N49	CO-N-N49		Neonate(<3kg)
LNOP Amtx	CO-T-A125	CO-F-A125	CO-N-A125		Adult(>30kg)
LNOP Pmtx	CO-T-P125	CO-F-P125	CO-N-P125		Pediatric(15-40kg)
LNOP Inf	CO-T-I125	CO-F-I125	CO-N-I125		Infant(1-15kg)
LNOP Neo	CO-T-N125	CO-F-N125	CO-N-N125		Neonate(<3kg)
M-LnCS Amtx	CO-T-A203	CO-F-A203	CO-N-A203		Adult(>30kg)
M-LNCS Pmtx	CO-T-P203	CO-F-P203	CO-N-P203		Pediatric(15-40kg)
M-LNCS Inf	CO-T-I203	CO-F-I203	CO-N-I203		Infant(1-15kg)
M-LnCS Neo	CO-T-N203	CO-F-N203	CO-N-N203		Neonate(<3kg)
<b>Nihon Kohden</b>					
TL-251T	CO-T-A16	CO-F-A16	CO-N-A16		Adult(>30kg)
TL-252T	CO-T-P16	CO-F-P16	CO-N-P16		Pediatric(15-40kg)
TL-252T	CO-T-I16	CO-F-I16	CO-N-I16		Infant(1-15kg)
TL-253T	CO-T-N16	CO-F-N16	CO-N-N16		Neonate(<3kg)

## 1.4 SPO2 Accessories

								
AC-01290	AC-01288	AC-01287	AC-01286	AC-01285	AC-01284	CO-1921	CO-1922	CO-1923
AC-01290	Ear sensor hanger (1 pcs)							
AC-01288	Adult ear clip Or animalsmall size lingua clip (1 pcs)							
AC-01287	Animal lingua clip,Big size (1 pcs)							
AC-01286	Animal lingua clip,Middle size (1 pcs)							
AC-01285	Reusable Silicon Y sensor wrap tape (1 pcs)							
AC-01284	Reusable 3 holes foam tape For All patient sizes (12 pcs Per bag)							
CO-1921	Disposable 3M adhesive textile tape For Adult and neonate (25 pcs Per bag)							
CO-1922	Disposable 3M adhesive textile tape for Pediatric and infant (25 pcs Per bag)							
CO-1923	Velcro foam wrap tape (12 pcs per bag)							